



Npy2r Antibody

Product Code	CSB-PA016036ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97295
Immunogen	Recombinant Mouse Neuropeptide Y receptor type 2 (Npy2r), partial (1-51AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide Y receptor type 2 (NPY2-R) (NPY-Y2 receptor) (Y2 receptor), Npy2r
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Gene Names	Npy2r
Product Modify	Unmodified