



AO440_004758 Antibody

Product Code	CSB-PA27269ZA01GJK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A0W0CTK3
Immunogen	Recombinant Candida glabrata AO440_004758 protein
Raised In	Rabbit
Species Reactivity	Candida glabrata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetyl-CoA carboxylase (AO440_004758)
Immunogen Species	Candida glabrata (Yeast) (Torulopsis glabrata)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	AO440_004758
Product Modify	Unmodified