



ANPEP Antibody

Product Code	CSB-PA001827ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15145
Immunogen	Recombinant Sus scrofa ANPEP protein
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aminopeptidase N (AP-N) (pAPN) (EC 3.4.11.2) (Alanyl aminopeptidase) (Aminopeptidase M) (AP-M) (Microsomal aminopeptidase) (gp130) (CD antigen CD13), ANPEP
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ANPEP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.