



ANPEP Antibody

Product Code	CSB-PA001827ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15145
Immunogen	Recombinant Pig Aminopeptidase N (ANPEP), partial (616-963AA)
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Aminopeptidase N (AP-N) (pAPN) (EC 3.4.11.2) (Alanyl aminopeptidase) (Aminopeptidase M) (AP-M) (Microsomal aminopeptidase) (gp130) (CD antigen CD13), ANPEP
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ANPEP
Product Modify	Unmodified