



MPN_372 Antibody

Product Code	CSB-PA300476ZA01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P75409
Immunogen	Recombinant Mycoplasma pneumoniae MPN_372 protein
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylating toxin CARDS (EC 2.4.2.-) (ADP-ribosyltransferase CARDS) (CARDX TX)
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MPN_372
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.