



BPHL Antibody

Product Code	CSB-PA774821ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WA6
Immunogen	Recombinant Human Valacyclovir hydrolase (BPHL) (38-291AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Valacyclovir hydrolase (VACVase) (Valacyclovirase) (EC 3.1.-.-) (Biphenyl hydrolase-like protein) (Biphenyl hydrolase-related protein) (Bph-rp) (Breast epithelial mucin-associated antigen) (MCNAA), BPHL, MCNAA
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Target Names	BPHL
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.