



sboA Antibody

Product Code	CSB-PA522484ZA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O07623
Immunogen	Recombinant Bacillus subtilis Subtilisin-A (sboA) (9-43AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Subtilisin-A (Antilisterial bacteriocin subtilisin), sboA, sbo
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order
Gene Names	sboA
Product Modify	Unmodified