



fpaP Antibody

Product Code	CSB-PA516921ZA01FEW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O05420
Immunogen	Recombinant Elizabethkingia meningoseptica fpaP protein
Raised In	Rabbit
Species Reactivity	Elizabethkingia meningoseptica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proline iminopeptidase (PIP) (EC 3.4.11.5) (Prolyl aminopeptidase) (PAP), fpaP
Immunogen Species	Elizabethkingia meningoseptica (Chryseobacterium meningosepticum)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fpaP
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.