




epiA Antibody

Product Code	CSB-PA357305ZA01FLL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08136
Immunogen	Recombinant Staphylococcus epidermidis Lantibiotic epidermin (epiA) (31-52AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus epidermidis
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Lantibiotic epidermin, epiA
Immunogen Species	Staphylococcus epidermidis
Target Names	epiA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: epiA antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 17.7 kDa

Observed band size: 20 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.