



nifH1 Antibody

Product Code	CSB-PA341970ZA01MEU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51602
Immunogen	Recombinant Methanobacterium ivanovii nifH1 protein
Raised In	Rabbit
Species Reactivity	Methanobacterium ivanovii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitrogenase iron protein 1 (EC 1.18.6.1) (Nitrogenase Fe protein 1) (Nitrogenase component II) (Nitrogenase reductase), nifH1
Immunogen Species	Methanobacterium ivanovii
Lead Time	Made-to-order (12-14 weeks)
Target Names	nifH1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.