



KMP-11 Antibody

Product Code	CSB-PA880045ZA01TQY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9U6Z1
Immunogen	Recombinant Trypanosoma cruzi KMP-11 protein
Raised In	Rabbit
Species Reactivity	Trypanosoma cruzi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetoplastid membrane protein 11 (KMP-11), KMP-11, KMP11
Immunogen Species	Trypanosoma cruzi
Lead Time	Made-to-order (12-14 weeks)
Target Names	KMP-11
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.