



Gucy2c Antibody

Product Code	CSB-PA010053ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23897
Immunogen	Recombinant Rattus norvegicus Gucy2c protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-stable enterotoxin receptor (STA receptor) (EC 4.6.1.2) (Guanylyl cyclase C) (GC-C) (Intestinal guanylate cyclase), Gucy2c, Guc2c
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Gucy2c
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.