



Susd4 Antibody

Product Code	CSB-PA805887ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BH32
Immunogen	Recombinant Mouse Sushi domain-containing protein 4 (Susd4), partial (42-319AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sushi domain-containing protein 4, Susd4
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order
Gene Names	Susd4
Product Modify	Unmodified