



ISM1 Antibody

Product Code	CSB-PA011850ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1AKI9
Immunogen	Recombinant Human Isthmin-1 (ISM1) (30-464AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Isthmin-1, ISM1, C20orf82 ISM
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Target Names	ISM1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.