




aur Antibody

Product Code	CSB-PA305800ZA01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81177
Immunogen	Recombinant Staphylococcus aureus Zinc metalloproteinase aureolysin (aur) (210-509AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc metalloproteinase aureolysin (EC 3.4.24.29) (Staphylococcus aureus neutral proteinase), aur
Immunogen Species	Staphylococcus aureus
Target Names	aur

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: aur antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39.3 kDa

Observed band size: 50 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.