



ctfB Antibody

Product Code	CSB-PA332931ZA01DUB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23673
Immunogen	Recombinant Clostridium acetobutylicum ctfB protein
Raised In	Rabbit
Species Reactivity	Clostridium acetobutylicum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Butyrate--acetoacetate CoA-transferase subunit B (Coat B) (EC 2.8.3.9) (Acetoacetyl-CoA:acetate/butyrate CoA-transferase subunit B), ctfB
Immunogen Species	Clostridium acetobutylicum (strain ATCC 824 / DSM 792 / JCM 1419 / LMG 5710 / VKM B-1787)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ctfB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.