



GPC Antibody

Product Code	CSB-PA362480ZA01LNP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08669
Immunogen	Recombinant Lassa virus GPC protein
Raised In	Rabbit
Species Reactivity	Lassa virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-glycoprotein polyprotein GP complex (Pre-GP-C) [Cleaved into: Stable signal peptide (SSP); Glycoprotein G1 (GP1); Glycoprotein G2 (GP2)], GPC, GP-C
Immunogen Species	Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	GPC
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.