



Tnmd Antibody

Product Code	CSB-PA024007ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9EP64
Immunogen	Recombinant Mus musculus Tnmd protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tenomodulin (TeM) (mTeM) (Chondromodulin-1-like protein) (ChM1L) (mChM1L) (Chondromodulin-I-like protein) (Myodulin) (Tendin), Tnmd, Chm1l
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Tnmd
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.