



Lims1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA859129ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99JW4 |
| Immunogen | Recombinant Mus musculus Lims1 protein |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | LIM and senescent cell antigen-like-containing domain protein 1 (Particularly interesting new Cys-His protein 1) (PINCH-1), Lims1, Pinch1 |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | Lims1 |
| Product Modify | Unmodified |