




# A27L Antibody

<b>Product Code</b>	CSB-PA26789ZA01VAR
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q89694
<b>Immunogen</b>	Recombinant Variola virus A27L protein (1-194aa )
<b>Species Reactivity</b>	Variola virus
<b>Tested Applications</b>	ELISA, WB
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Variola virus (isolate Human/India/Ind3/1967) (VARV) (Smallpox virus)
<b>Research Area</b>	Others
<b>Target Names</b>	A27L

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: A27L antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.