




A27L Antibody

Product Code	CSB-PA26789ZA01VAR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q89694
Immunogen	Recombinant Variola virus A27L protein (1-194aa)
Species Reactivity	Variola virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Variola virus (isolate Human/India/Ind3/1967) (VARV) (Smallpox virus)
Research Area	Others
Target Names	A27L

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: A27L antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa