



Os09g0401200 Antibody

Product Code	CSB-PA732309ZA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ES52
Immunogen	Recombinant Oryza sativa subsp. japonica Os09g0401200 protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	TPR repeat-containing thioredoxin TDX (OsTrx26) (Tetratricoredoxin) (OsTDX)
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Os09g0401200
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.