



tmpB Antibody

Product Code	CSB-PA333469ZA01THS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29720
Immunogen	Recombinant Treponema phagedenis tmpB protein
Raised In	Rabbit
Species Reactivity	Treponema phagedenis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Treponemal membrane protein B (Antigen TmpB), tmpB
Immunogen Species	Treponema phagedenis
Lead Time	Made-to-order (12-14 weeks)
Gene Names	tmpB
Product Modify	Unmodified