



uvvE Antibody

Product Code	CSB-PA467668ZA01BQP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7IR00
Immunogen	Recombinant Bacillus cereus uvvE protein
Raised In	Rabbit
Species Reactivity	Bacillus cereus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UV DNA damage endonuclease (UV-endonuclease) (UVED) (EC 3.-.-.-), uvvE
Immunogen Species	Bacillus cereus (strain G9842)
Lead Time	Made-to-order (12-14 weeks)
Target Names	uvvE
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.