



troA Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA309098ZA01TPR |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P96116 |
| Immunogen | Recombinant Treponema pallidum troA protein |
| Raised In | Rabbit |
| Species Reactivity | Treponema pallidum |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Periplasmic zinc-binding protein TroA (Tromp-1), troA, troMP1 |
| Immunogen Species | Treponema pallidum (strain Nichols) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | troA |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |