



HSP100 Antibody

Product Code	CSB-PA517681ZA01TQY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15885
Immunogen	Recombinant Trypanosoma cruzi HSP100 protein
Raised In	Rabbit
Species Reactivity	Trypanosoma cruzi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein 100 (Protein CLP) (Fragment), HSP100
Immunogen Species	Trypanosoma cruzi
Lead Time	Made-to-order (12-14 weeks)
Target Names	HSP100
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.