



UBC5A Antibody

Product Code	CSB-PA848293ZA01OFG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8S920
Immunogen	Recombinant Oryza sativa subsp. japonica UBC5A protein
Raised In	Rabbit
Species Reactivity	Oryza sativa subsp. japonica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 5A (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme 5A) (Ubiquitin carrier protein 5a) (OsUBC5a) (Ubiquitin-protein ligase 5A), UBC5A
Immunogen Species	Oryza sativa subsp. japonica (Rice)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UBC5A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.