



GP Antibody

Product Code	CSB-PA330440ZA01LBR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35254
Immunogen	Recombinant Lake Victoria marburgvirus GP protein
Raised In	Rabbit
Species Reactivity	Lake Victoria marburgvirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein (GP1,2) (GP) (Virion spike glycoprotein) [Cleaved into: GP1; GP2], GP
Immunogen Species	Lake Victoria marburgvirus (strain Popp-67) (MARV) (Marburg virus (strain West Germany/Popp/1967))
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GP
Product Modify	Unmodified