



BQ2027_MB3645C Antibody

Product Code	CSB-PA353261ZA01MVH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P65088
Immunogen	Recombinant Mycobacterium bovis BQ2027_MB3645C protein
Raised In	Rabbit
Species Reactivity	Mycobacterium bovis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	EspC protein homolog
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Lead Time	Made-to-order (12-14 weeks)
Target Names	BQ2027_MB3645C
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.