



tpd Antibody

Product Code	CSB-PA322536ZA01TPR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19478
Immunogen	Recombinant Treponema pallidum tpd protein
Raised In	Rabbit
Species Reactivity	Treponema pallidum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	34 kDa membrane antigen (Pathogen-specific membrane antigen), tpd
Immunogen Species	Treponema pallidum (strain Nichols)
Lead Time	Made-to-order (12-14 weeks)
Target Names	tpd
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.