



hlgC Antibody

Product Code	CSB-PA751944ZA01SKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7A3S2
Immunogen	Recombinant Staphylococcus aureus Gamma-hemolysin component C (hlgC) (30-315AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-hemolysin component C, hlgC
Immunogen Species	Staphylococcus aureus (strain N315)
Lead Time	Made-to-order
Gene Names	hlgC
Product Modify	Unmodified