



PAPD5 Antibody

Product Code	CSB-PA843329ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDF8
Immunogen	Recombinant Human Non-canonical poly (A) RNA polymerase PAPD5 (PAPD5) (1-572AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-canonical poly(A) RNA polymerase PAPD5 (EC 2.7.7.19) (PAP-associated domain-containing protein 5) (Terminal uridylyltransferase 3) (TUTase 3) (Topoisomerase-related function protein 4-2) (TRF4-2), PAPD5
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	PAPD5
Product Modify	Unmodified