



sar Antibody

Product Code	CSB-PA365461ZA01APN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00655
Immunogen	Recombinant Aspergillus giganteus sar protein
Raised In	Rabbit
Species Reactivity	Aspergillus giganteus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease alpha-sarcin (EC 3.1.27.10) (rRNA endonuclease), sar
Immunogen Species	Aspergillus giganteus
Lead Time	Made-to-order (12-14 weeks)
Target Names	sar
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.