



Lpl Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA013065ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P11152 |
| Immunogen | Recombinant Mus musculus Lpl protein |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lipoprotein lipase (LPL) (EC 3.1.1.34), Lpl |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | Lpl |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |