



alp1 Antibody

Product Code	CSB-PA329798ZA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28296
Immunogen	Recombinant Neosartorya fumigata alp1 protein
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkaline protease 1 (ALP) (EC 3.4.21.63) (Aspergillopeptidase B) (Aspergillus proteinase B) (Elastase) (Elastinolytic serine proteinase) (Oryzin) (allergen Asp f 13), alp1, alk1 alp
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	alp1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.