



# Adamts13 Antibody

<b>Product Code</b>	CSB-PA745079ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q769J6
<b>Immunogen</b>	Recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 13 (Adamts13), partial (904-1137AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein A/G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	A disintegrin and metalloproteinase with thrombospondin motifs 13 (ADAM-TS 13) (ADAM-TS13) (ADAMTS-13) (EC 3.4.24.87) (von Willebrand factor-cleaving protease) (vWF-CP) (vWF-cleaving protease), Adamts13, Gm710
<b>Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order
<b>Gene Names</b>	Adamts13
<b>Product Modify</b>	Unmodified