



# Nr2f6 Antibody

<b>Product Code</b>	CSB-PA016058ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P43136
<b>Immunogen</b>	Recombinant Mus musculus Nr2f6 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Nuclear receptor subfamily 2 group F member 6 (COUP transcription factor 3) (COUP-TF3) (V-erbA-related protein 2) (EAR-2), Nr2f6, Ear-2 Ear2 Erbal2
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	Nr2f6
<b>Product Modify</b>	Unmodified