



Exochitosanase Antibody

Product Code	CSB-PA306651ZA01AVD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83488
Immunogen	Recombinant Aspergillus flavus Exochitosanase protein
Raised In	Rabbit
Species Reactivity	Aspergillus flavus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exochitosanase (EC 3.2.1.-) (Exo-beta-D-glucosaminidase) (Fragment)
Immunogen Species	Aspergillus flavus
Lead Time	Made-to-order (12-14 weeks)
Target Names	Exochitosanase
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.