



Serine protease inhibitor Antibody

Product Code	CSB-PA305524ZA01LDV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81639
Immunogen	Recombinant Lentinula edodes Serine protease inhibitor protein
Raised In	Rabbit
Species Reactivity	Lentinula edodes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor
Immunogen Species	Lentinula edodes (Shiitake mushroom) (Lentinus edodes)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.