



dsbA Antibody

Product Code	CSB-PA321092ZA01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13320
Immunogen	Recombinant Enterobacteria phage T4 dsbA protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Double-stranded DNA-binding protein (DsDNA-binding protein A), dsbA, rpbB
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dsbA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.