



Ctsl Antibody

Product Code	CSB-PA006193ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06797
Immunogen	Recombinant Mus musculus Ctsl protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin L1 (EC 3.4.22.15) (Cathepsin L) (Major excreted protein) (MEP) (p39 cysteine proteinase) [Cleaved into: Cathepsin L1 heavy chain; Cathepsin L1 light chain], Ctsl, Ctsl1
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Ctsl
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.