



MMP20 Antibody

Product Code	CSB-PA014667ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60882
Immunogen	Recombinant Homo sapiens MMP20 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix metalloproteinase-20 (MMP-20) (EC 3.4.24.-) (Enamel metalloproteinase) (Enamelysin), MMP20
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	MMP20
Product Modify	Unmodified