



btuD Antibody

Product Code	CSB-PA538512ZA01HTL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0R5G4
Immunogen	Recombinant Halobacterium salinarum btuD protein
Raised In	Rabbit
Species Reactivity	Halobacterium salinarum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cobalamin import ATP-binding protein BtuD (EC 3.6.3.33) (Vitamin B12-transporting ATPase), btuD
Immunogen Species	Halobacterium salinarum (strain ATCC 29341 / DSM 671 / R1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	btuD
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.