



acpS Antibody

Product Code	CSB-PA669769ZA01SBAF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q48RM7
Immunogen	Recombinant Streptococcus pyogenes serotype M28 acpS protein
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes serotype M28
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Holo-[acyl-carrier-protein] synthase (Holo-ACP synthase) (EC 2.7.8.7) (4'-phosphopantetheinyl transferase AcpS), acpS
Immunogen Species	Streptococcus pyogenes serotype M28 (strain MGAS6180)
Lead Time	Made-to-order (12-14 weeks)
Target Names	acpS
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.