



4 Antibody

Product Code	CSB-PA356120ZA01EEB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03692
Immunogen	Recombinant Enterobacteria phage T7 4 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T7
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA primase/helicase (EC 2.7.7.-) (EC 3.6.4.12) (Gene product 4) (Gp4)
Immunogen Species	Enterobacteria phage T7 (Bacteriophage T7)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	4
Product Modify	Unmodified