



METR N Antibody

Product Code	CSB-PA892322ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJH8
Immunogen	Recombinant Human Meteorin (METRN) (24-293AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Meteorin, METRN, C16orf23
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	METRN
Product Modify	Unmodified