




DBP Antibody

Product Code	CSB-PA365892ZA01HIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03265
Immunogen	Recombinant Human adenovirus C serotype 5 DNA-binding protein (DBP), partial (174-529AA)
Raised In	Rabbit
Species Reactivity	Human adenovirus C serotype 5
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein (DBP) (Early 2A protein) (Early E2A DNA-binding protein), DBP
Immunogen Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Target Names	DBP

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: DBP antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43.5 kDa

Observed band size: 49 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.