



btuD Antibody

Product Code	CSB-PA538512ZA01HTL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0R5G4
Immunogen	Recombinant Halobacterium salinarum Cobalamin import ATP-binding protein BtuD (btuD) (1-398AA)
Raised In	Rabbit
Species Reactivity	Halobacterium salinarum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cobalamin import ATP-binding protein BtuD (EC 3.6.3.33) (Vitamin B12-transporting ATPase), btuD
Immunogen Species	Halobacterium salinarum (strain ATCC 29341 / DSM 671 / R1)
Lead Time	Made-to-order
Gene Names	btuD
Product Modify	Unmodified