



# DPYSL3 Antibody

<b>Product Code</b>	CSB-PA622756ZA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q14195
<b>Immunogen</b>	Recombinant Homo sapiens DPYSL3 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Homo sapiens
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Dihydropyrimidinase-related protein 3 (DRP-3) (Collapsin response mediator protein 4) (CRMP-4) (Unc-33-like phosphoprotein 1) (ULIP-1), DPYSL3, CRMP4, DRP3, ULIP, ULIP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	DPYSL3
<b>Product Modify</b>	Unmodified