



# crf1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA818279ZA01NGS   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8J0P4  |
| <b>Immunogen</b>           | Recombinant Neosartorya fumigata crf1 protein   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Neosartorya fumigata  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                      |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Probable glycosidase crf1 (EC 3.2.-.-) (Crh-like protein 1) (allergen Asp f 9), crf1                  |
| <b>Immunogen Species</b>   | Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus) |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)   |
| <b>Target Names</b>        | crf1  |
| <b>Product Modify</b>      | Unmodified  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                           |